

Jones

Medical Supply

Home Oxygen & Medical Equipment

Thank you for allowing Jones Medical Supply to provide your home oxygen equipment. Below you will find step by step instructions on how to use your Invacare oxygen concentrator.

Important reminder; no smoking or open flames within 25 feet of you or the concentrator.



Equipment may vary, but in most cases your equipment will look like one of these two Invacare products.



- After plugging your machine into a safe wall outlet, free from any curtains or loose bedding (that can be pulled into the filters), press the black power button. **DO NOT** use surge protectors or extension cords.
- We've provided a tubing and cannulas for your concentrator.
- Connect your tubing to the brass or black port, ensure it's secure.
- Adjust the liter flow to the prescribed rate (which is indicated on your paperwork) by turning the dial and watching the ball move in the cylinder.
- Make sure there are no obstacles between you and the machine, **NO OPEN FLAMES** of any type, no rugs, furniture or other potential tripping hazards.

